



## Stacked Large Capacity Shaking Incubator

### Product Introduction


The BJPX-SDT series stackable tabletop thermostatic shakers are widely used in scenarios requiring strict control over temperature and oscillation frequency, such as bacterial cultivation, fermentation, hybridization, biochemical reactions, and research on enzymes and cell tissues. They have wide and critical applications in research fields including medicine, biology, molecular biology, pharmaceuticals, food science, and environmental protection.

These shakers are suitable for cell culture, extraction experiments, sample mixing, staining and destaining, hybridization experiments, solubility tests, plasmid purification, bacterial suspension preparation, and elution processes.

### Fixture Display



Universal spring clip



Test Tube Rack



Rod-Clamp Shaking Plate

## Application Scenarios



Extraction



Pharmacy



Hospital

## Product Features



PID microcomputer intelligent temperature controller with precise temperature control.



Electrically heated tempered glass observation window with anti-condensation function



5.0-inch touchscreen for simple and intuitive parameter setting



Equipped with UV sterilization function to keep the chamber clean and sterile



Standard with LED lighting function, which automatically turns on when the chamber door is opened



It supports up to three stackable units for space-saving design, features a convenient pull-out tray, and enables flexible operation.



High-quality DC brushless motor with low vibration and low noise



Rotating speed adjustable from 30 to 300 (r/min)



Imported compressor, air-cooled R290 non-CFC refrigerant



Temperature control range: RT+5°C ~60°C / 4°C~60°C



Curve programming settings: cyclic/stepwise/eight-segment curve programming



Overtemperature alarm, abnormal shutdown alarm, motor fault alarm, door open timeout alarm

## Fixture Selection

Product Model	Clamp Model							
BJPX-SDW10 BJPX-SDW20	Maximum capacity of Erlenmeyer flask					Test tube rack	96-well plate	Rod clamp shaking plate
	250ml	500ml	1000ml	2000ml	3000ml			
	34个	24个	15个	8个	6个			