



Water Bath Incubator Shaker

Product Introduction

Water Bath Incubator Shaker

Combines temperature-controlled water bath with orbital shaking for precision culturing.

Essential Equipment

For plant/biological/microbiological research, genetic/viral studies, medical/environmental sciences, and production sectors.

Product Features



Premium Digital Controller
Drift-free speed stability



Splash-Proof Baffle, Eliminates water splashing at high speeds



Reciprocating Motion
22mm Amplitude



Rapid Temp/Speed Adjustment
Saves valuable experiment time



Powder-Coated Steel Housing
Durable and rugged construction



PID Microprocessor Precision Thermal Control



Programmable Timer
Multiple timing modes
Spring-Loaded Clamps



Universal compatibility with flasks, beakers & irregular samples



Over-Temp & Over-Speed Alarms



Parameter Memory
Auto-resume after power recovery



Temp Range: RT+5° C to 100° C



Adjustable Speed: 30-300 rpm