

Water Bath



Product Introduction

Water Bath
Essential for laboratories across universities, research institutes, and corporate R&D facilities.
Ideal for:
Sample distillation • Drying/concentration • Soaking procedures • Serum testing • Biochemical incubations • Sterilization of instruments via boiling

Product Features

- Cold-Rolled Steel Shell
Enhanced corrosion resistance
 - Premium Stainless Steel Inner Tank
Single-stamped construction
 - LED Display
Real-time temp & countdown
 - Retractable Drain
Sleek, durable & convenient
- 5-RING SPECS (MM):
Φ37 • Φ58 • Φ79 • Φ100 • Φ120
 - Built-in Temp Calibration
 - Auto Dry-Burn Protection
Timer: 1-999 min
 - Parameter Memory
Auto-resume after power recovery

Product selection

Models	HH-S1	HH-S2	HH-S4	HH-S6	HH-S8
Volume (L)	3.5	6.8	13	20	26
Temp Range (°C)	RT +5-99.9				
Comments	Single hole	Single column two holes	Single column four holes	Single column six holes	Single column eight holes