













## Vortex Mixer

### Product Introduction

The bjpx-VW Vortex Mixer rapidly homogenizes liquids, powders, and other samples using high-speed vortex motion.

This versatile instrument is widely employed in laboratories and testing facilities across multiple sectors, including environmental monitoring, healthcare, petrochemicals, food processing, metallurgy, PCR preparation areas, as well as in academic research and industrial production settings.

### Product Features

-  Crafted with premium ABS, durability and refined aesthetics.
-  Brushless DC Motor: Seamless Speed Control, Quiet Operation, Extended Lifespan.
-  Digital display, precise control of time and speed.
-  Fast response and speed adjustments.
-  Robust construction, superior vibration damping for high-speed operation.
-  Multiple optional models and kits
-  Crafted with premium ABS, durability and refined aesthetics.
-  Adjustable speed: 20-3000 rpm
-  Orbital Oscillation: 4mm Amplitude
-  Pressure-sensing control with pulse/-continuous modes